

SemenAssay® POX -Semen Leukocytes Test Kit (Peroxidase Staining)

The kit is used for quantitative detection of leukocyte concentration in human semen specimen and screening the male reproductive tract infection.

The sketch map of test results:



The cells that were stained brown are peroxidase-positive.



There are 40 vials of reaction A per kit and the reaction A is directly used for testing conveniently.

Features:

- The only product contains the reagents of two components and it has less operation procedures than the competitive products. It is convenient for clinical use.
- It is ready for use and easy to handle.
- The HRP control solution is used for quality control and guarantees precision of the results.
- 40 tests per kit and meet the need of batch testing.
- CE approved

