

SpermFunc® Diff-Quik - Staining Kit for Spermatozoan Morphology (Diff-Quik Rapid Staining Method)

The kit is mainly used for evaluating the characteristics of human sperm morphology. It can also be used for the detection of sperm morphology of other mammals.

Features:

- A mature method recommended by WHO.
- Ready-to-use and faster than conventional modified Papanicolaou stain method.
- Special slide barrel prepared for staining. Simple and easy-to-use.
- Shelf life up to 2 years.
- 3 x 100ml or 3 x 500ml per kit. More Choices for lab requirement. Meet the need of batch testing.
- CE approved

The sketch maps of test results by Diff-Quik:

